

Anti-CD94 Mouse Monoclonal Antibody, PE Conjugated



Catalog #: F0566-PE

Aliases

KLRD1; Killer Cell Lectin Like Receptor D1; CD94; Killer Cell Lectin-Like Receptor Subfamily D, Member 1; Natural Killer Cells Antigen CD94; NK Cell Receptor; CD94 Antigen; KP43; Killer Cell Lectin-Like Receptor Subfamily D Member 1

Background

Gene Name: KLRD1

NCBI Gene Entry: [3824](#)

UniProt Entry: [Q13241](#)

Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB14790

Species Reactivity: Human, african green, baboon, cynomolgus

Applications Tested: Flow Cytometry (FCM)

Conjugation: PE

Immunogen

The immunogen information is proprietary.

Isotype

Mouse IgG1, κ

Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

Storage

Store at 4 °C and protected from prolonged exposure to light.

Recommended Dilutions

5 μ L/Tests

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041


ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-CD94 Mouse Monoclonal Antibody, PE Conjugated

PAGE 2



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM