

# Anti-CD95 Mouse Monoclonal Antibody, PE Conjugated



**Catalog #: F0567-PE**

## Aliases

FAS; Fas Cell Surface Death Receptor; TNFRSF6; APO-1; CD95; APT1; FAS1; Tumor Necrosis Factor Receptor Superfamily, Member 6; Tumor Necrosis Factor Receptor Superfamily Member 6; Fas (TNF Receptor Superfamily, Member 6); Apoptosis-Mediating Surface Antigen FAS; TNF Receptor Superfamily Member 6; FASLG Receptor; CD95 Antigen; Mutant Tumor Necrosis Receptor Superfamily Member 6; Apoptosis Signaling Receptor FAS; APO-1 Cell Surface Antigen; Apoptosis Antigen 1; Apo-1 Antigen; Fas AMA; ALPS1A; FASTM

## Background

Gene Name: FAS

NCBI Gene Entry: [355](#)

UniProt Entry: [P25445](#)

## Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB14795

Species Reactivity: Human, african green, baboon, capuchin monkey, chimpanzee, marmoset, cynomolgus, rhesus

Applications Tested: Flow Cytometry (FCM)

Conjugation: PE

## Immunogen

The immunogen information is proprietary.

## Isotype

Mouse IgG1,  $\kappa$

## Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

## Storage

Store at 4 °C and protected from prolonged exposure to light.

## Recommended Dilutions

5  $\mu$ L/Tests

---

### SUPPORT

SUPPORT@GENUINBIOTECH.COM  
TEL: +1-540-855-7041

### ORDERS

SALES@GENUINBIOTECH.COM  
FAX: +1-540-855-7041


[WWW.GENUINBIOTECH.COM](http://WWW.GENUINBIOTECH.COM)

# Anti-CD95 Mouse Monoclonal Antibody, PE Conjugated

PAGE 2

**Note:** This product is for research use only.

## Validation Data



### SUPPORT

SUPPORT@GENUINBIOTECH.COM  
TEL: +1-540-855-7041

### ORDERS

SALES@GENUINBIOTECH.COM  
FAX: +1-540-855-7041

[WWW.GENUINBIOTECH.COM](http://WWW.GENUINBIOTECH.COM)