

Anti-CD109 Mouse Monoclonal Antibody, PE Conjugated



Catalog #: F0575-PE

Aliases

CD109; CD109 Molecule; CPAMD7; C3 And PZP-Like Alpha-2-Macroglobulin Domain-Containing Protein 7; 150 KDa TGF-Beta-1-Binding Protein; Platelet-Specific Gov Antigen; DKFZp762L1111; CD109 Antigen; FLJ38569; P180; R150; CD109 Antigen (Gov Platelet Alloantigens); Activated T-Cell Marker CD109; Gov Platelet Alloantigens

Background

Gene Name: CD109

NCBI Gene Entry: [135228](#)

UniProt Entry: [Q6YHK3](#)

Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB15475

Species Reactivity: Human

Applications Tested: Flow Cytometry (FCM)

Conjugation: PE

Immunogen

The immunogen information is proprietary.

Isotype

Mouse IgG1

Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

Storage

Store at 4 °C and protected from prolonged exposure to light.

Recommended Dilutions

5 µL/Tests

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

