

Anti-CD133 Mouse Monoclonal Antibody, PE Conjugated



Catalog #: F0581-PE

Aliases

PROM1; Prominin 1; CORD12; PROML1; AC133; CD133; RP41; Retinitis Pigmentosa 41, Cone-Rod Dystrophy 12; Prominin-Like Protein 1; AC133 Antigen; Antigen AC133; Prominin-1; MCDR2; STGD4; Stargardt Disease 4 (Autosomal Dominant); Hematopoietic Stem Cell Antigen; Macular Dystrophy, Retinal 2; Prominin (Mouse)-Like 1; CD133 Antigen; HProminin; MSTP061

Background

Gene Name: PROM1

NCBI Gene Entry: [8842](#)

UniProt Entry: [O43490](#)

Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB14845

Species Reactivity: Human

Applications Tested: Flow Cytometry (FCM)

Conjugation: PE

Immunogen

The immunogen information is proprietary.

Isotype

Mouse IgG1, κ

Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

Storage

Store at 4 °C and protected from prolonged exposure to light.

Recommended Dilutions

5 μ L/Tests

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041


ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-CD133 Mouse Monoclonal Antibody, PE Conjugated

PAGE 2



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM