

Anti-CD140a Mouse Monoclonal Antibody, PE Conjugated



Catalog #: F0583-PE

Aliases

PDGFRA; Platelet Derived Growth Factor Receptor Alpha; PDGFR2; CD140a; Platelet-Derived Growth Factor Receptor, Alpha Polypeptide; Alpha-Type Platelet-Derived Growth Factor Receptor; Platelet-Derived Growth Factor Receptor Alpha; Platelet-Derived Growth Factor Receptor 2; CD140 Antigen-Like Family Member A; CD140a Antigen; PDGF-R-Alpha; EC 2.7.10.1; PDGFR-2; GAS9; Alpha Platelet-Derived Growth Factor Receptor; Platelet-Derived Growth Factor Alpha Receptor; PDGFR-Alpha; RHEPDGFRA; EC 2.7.10

Background

Gene Name: PDGFRA

NCBI Gene Entry: [5156](#)

UniProt Entry: [P16234](#)

Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB14860

Species Reactivity: Human

Applications Tested: Flow Cytometry (FCM)

Conjugation: PE

Immunogen

The immunogen information is proprietary.

Isotype

Mouse IgG1,k

Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

Storage

Store at 4 °C and protected from prolonged exposure to light.

Recommended Dilutions

5 µL/Tests

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041


WWW.GENUINBIOTECH.COM

Anti-CD140a Mouse Monoclonal Antibody, PE Conjugated

PAGE 2

Note: This product is for research use only.

Validation Data



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM