

# Anti-CD140b Mouse Monoclonal Antibody, PE Conjugated



**Catalog #: F0584-PE**

## Aliases

PDGFRB; Platelet Derived Growth Factor Receptor Beta; PDGFR1; CD140b; JTK12; PDGFR; Platelet-Derived Growth Factor Receptor, Beta Polypeptide; Beta-Type Platelet-Derived Growth Factor Receptor; Platelet-Derived Growth Factor Receptor Beta; Platelet-Derived Growth Factor Receptor 1; CD140 Antigen-Like Family Member B; PDGF-R-Beta; EC 2.7.10.1; PDGFR-Beta; PDGFR-1; Beta Platelet-Derived Growth Factor Receptor; Activated Tyrosine Kinase PDGFRB; CD140b Antigen; NDEL1-PDGFRB; EC 2.7.10; IBGC4; PENTT; IMF1; KOGS

## Background

Gene Name: PDGFRB

NCBI Gene Entry: [5159](#)

UniProt Entry: [P09619](#)

## Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB14865

Species Reactivity: Human

Applications Tested: Flow Cytometry (FCM)

Conjugation: PE

## Immunogen

The immunogen information is proprietary.

## Isotype

Mouse IgG1,k

## Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

## Storage

Store at 4 °C and protected from prolonged exposure to light.

## Recommended Dilutions

5 µL/Tests

---

### SUPPORT

SUPPORT@GENUINBIOTECH.COM  
TEL: +1-540-855-7041

### ORDERS

SALES@GENUINBIOTECH.COM  
FAX: +1-540-855-7041


[WWW.GENUINBIOTECH.COM](http://WWW.GENUINBIOTECH.COM)

# Anti-CD140b Mouse Monoclonal Antibody, PE Conjugated

PAGE 2

**Note:** This product is for research use only.

## Validation Data



### SUPPORT

SUPPORT@GENUINBIOTECH.COM  
TEL: +1-540-855-7041

### ORDERS

SALES@GENUINBIOTECH.COM  
FAX: +1-540-855-7041

[WWW.GENUINBIOTECH.COM](http://WWW.GENUINBIOTECH.COM)