

Anti-CD184 Mouse Monoclonal Antibody, PE Conjugated



Catalog #: F0592-PE

Aliases

CXCR4; C-X-C Motif Chemokine Receptor 4; LESTR; Fusin; HM89; D2S201E; NPY3R; HSY3RR; NPY3R; CD184; NPYR; Leukocyte-Derived Seven Transmembrane Domain Receptor; Lipopolysaccharide-Associated Protein 3; Stromal Cell-Derived Factor 1 Receptor; Chemokine (C-X-C Motif) Receptor 4; C-X-C Chemokine Receptor Type 4; LPS-Associated Protein 3; SDF-1 Receptor; CD184 Antigen; LAP-3; NPYRL; FB22; LCR1; Seven-Transmembrane-Segment Receptor, Spleen; Chemokine (C-X-C Motif), Receptor 4 (Fusin); Seven Transmembrane Helix Receptor; Neuropeptide Y Receptor Y3; Neuropeptide Y3 Receptor; WHIMS1; CXC-R4; CXCR-4; WHIMS; LAP3; WHIM

Background

Gene Name: CXCR4

NCBI Gene Entry: [7852](#)

UniProt Entry: [P61073](#)

Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB14900

Species Reactivity: Human

Applications Tested: Flow Cytometry (FCM)

Conjugation: PE

Immunogen

The immunogen information is proprietary.

Isotype

Mouse IgG2a, κ

Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

Storage

Store at 4 °C and protected from prolonged exposure to light.

Recommended Dilutions

5 μ L/Tests

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041


WWW.GENUINBIOTECH.COM

Anti-CD184 Mouse Monoclonal Antibody, PE Conjugated

PAGE 2

Note: This product is for research use only.

Validation Data



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM