

# Anti-CD200 Mouse Monoclonal Antibody, PE Conjugated



**Catalog #: F0593-PE**

## Aliases

CD200; CD200 Molecule; CD200 Antigen; OX-2; MOX1; MOX2; MRC; Antigen Identified By Monoclonal Antibody MRC OX-2; OX-2 Membrane Glycoprotein

## Background

Gene Name: CD200

NCBI Gene Entry: [4345](#)

UniProt Entry: [P41217](#)

## Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB14905

Species Reactivity: Human, african green, baboon

Applications Tested: Flow Cytometry (FCM)

Conjugation: PE

## Immunogen

The immunogen information is proprietary.

## Isotype

Mouse IgG1,  $\kappa$

## Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

## Storage

Store at 4 °C and protected from prolonged exposure to light.

## Recommended Dilutions

5  $\mu$ L/Tests

**Note:** This product is for research use only.

## Validation Data

---

### SUPPORT

SUPPORT@GENUINBIOTECH.COM  
TEL: +1-540-855-7041


### ORDERS

SALES@GENUINBIOTECH.COM  
FAX: +1-540-855-7041

[WWW.GENUINBIOTECH.COM](http://WWW.GENUINBIOTECH.COM)

## Anti-CD200 Mouse Monoclonal Antibody, PE Conjugated

PAGE 2



### SUPPORT

SUPPORT@GENUINBIOTECH.COM  
TEL: +1-540-855-7041

### ORDERS

SALES@GENUINBIOTECH.COM  
FAX: +1-540-855-7041

[WWW.GENUINBIOTECH.COM](http://WWW.GENUINBIOTECH.COM)