

Anti-CD273 Mouse Monoclonal Antibody, PE Conjugated



Catalog #: F0597-PE

Aliases

PDCD1LG2; Programmed Cell Death 1 Ligand 2; PD-L2; CD273; B7-DC; PDL2; B7DC; BA574F11.2; Btdc; B7 Dendritic Cell Molecule; Programmed Death Ligand 2; Butyrophilin B7-DC; PDCD1 Ligand 2; PDCD1L2; PD-1-Ligand 2; PD-1 Ligand 2; CD273 Antigen

Background

Gene Name: PDCD1LG2

NCBI Gene Entry: [80380](#)

UniProt Entry: [Q9BQ51](#)

Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB14925

Species Reactivity: Human

Applications Tested: Flow Cytometry (FCM)

Conjugation: PE

Immunogen

The immunogen information is proprietary.

Isotype

Mouse IgG2a, κ

Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

Storage

Store at 4 °C and protected from prolonged exposure to light.

Recommended Dilutions

5 μ L/Tests

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041


ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-CD273 Mouse Monoclonal Antibody, PE Conjugated

PAGE 2



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM