

Anti-HLA-DQ Mouse Monoclonal Antibody, PE Conjugated



Catalog #: F0607-PE

Aliases

HLA-DQA1; Major Histocompatibility Complex, Class II, DQ Alpha 1; CELIAC1; HLA Class II Histocompatibility Antigen, DQ Alpha 1 Chain; MHC Class II DQA1; DC-1 Alpha Chain; DC-Alpha; HLA-DQA; HLA-DCA; HLA Class II Histocompatibility Antigen DQ Alpha Chain; MHC Class II DQ Alpha Chain; MHC Class II HLA-DQ-Alpha-1; MHC Class II Antigen DQA1; MHC Class II Protein; MHC HLA-DQ Alpha; HLA-DQA1*; HLA-DQB1; DQ-A1; DQA1; HLA-DQB1; Major Histocompatibility Complex, Class II, DQ Beta 1; CELIAC1; HLA-DQB; IDDM1; HLA Class II Histocompatibility Antigen, DQ Beta 1 Chain; MHC Class II Antigen DQB1; MHC Class II HLA-DQ Beta Glycoprotein; MHC Class II Antigen HLA-DQ-Beta-1; MHC Class II DQ Beta Chain

Background

Gene Name: HLA-DQA1/HLA-DQB1

NCBI Gene Entry: [3117/3119](#)

UniProt Entry: [P01909/P01920](#)

Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB15330

Species Reactivity: Human

Applications Tested: Flow Cytometry (FCM)

Conjugation: PE

Immunogen

The immunogen information is proprietary.

Isotype

Mouse IgG2a

Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

Storage

Store at 4 °C and protected from prolonged exposure to light.

Recommended Dilutions

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM


Anti-HLA-DQ Mouse Monoclonal Antibody, PE Conjugated

PAGE 2

5 µL/Tests

Note: This product is for research use only.

Validation Data



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM