

# Anti-FOXP3 Mouse Monoclonal Antibody, PE Conjugated



**Catalog #: F0609-PE**

## Aliases

FOXP3; Forkhead Box P3; DIETER; XPID; AIID; PIDX; IPEX; JM2; Immune Dysregulation, Polyendocrinopathy, Enteropathy, X-Linked; Forkhead Box Protein P3; SCURFIN; Immunodeficiency, Polyendocrinopathy, Enteropathy, X-Linked; FOXP3delta7; Scurfin

## Background

Gene Name: FOXP3

NCBI Gene Entry: [50943](#)

UniProt Entry: [Q9BZS1](#)

## Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB15490

Species Reactivity: Human

Applications Tested: Intracellular Staining for Flow Cytometry ☐ ICFC ☐

Conjugation: PE

## Immunogen

The immunogen information is proprietary.

## Isotype

Mouse IgG1,  $\kappa$

## Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

## Storage

Store at 4 °C and protected from prolonged exposure to light.

## Recommended Dilutions

5  $\mu$ L/Tests

**Note:** This product is for research use only.

## Validation Data

---

### SUPPORT

SUPPORT@GENUINBIOTECH.COM  
TEL: +1-540-855-7041

### ORDERS

SALES@GENUINBIOTECH.COM  
FAX: +1-540-855-7041

[WWW.GENUINBIOTECH.COM](http://WWW.GENUINBIOTECH.COM)

