

Anti-IL-4 Rat Monoclonal Antibody, PE Conjugated



Catalog #: F0615-PE

Aliases

IL4; Interleukin 4; IL-4; Lymphocyte Stimulatory Factor 1; BCGF-1; BCGF1; BSF1; B_cell Stimulatory Factor 1; B Cell Growth Factor 1; Interleukin-4; Binetrakin; Pittrakinra; MGC79402; BSF-1; B-Cell Stimulatory Factor 1

Background

Gene Name: IL4

NCBI Gene Entry: [3565](#)

UniProt Entry: [P05112](#)

Application Information

Clonality: Rat monoclonal antibody

Clone ID: 24GB14985

Species Reactivity: Human

Applications Tested: Intracellular Staining for Flow Cytometry ☐ ICFC ☐

Conjugation: PE

Immunogen

The immunogen information is proprietary.

Isotype

Rat IgG1, κ

Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

Storage

Store at 4 °C and protected from prolonged exposure to light.

Recommended Dilutions

5 μ L/Tests

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041


ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-IL-4 Rat Monoclonal Antibody, PE Conjugated

PAGE 2



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM