

Anti-IL-6 Rat Monoclonal Antibody, PE Conjugated



Catalog #: F0616-PE

Aliases

IL6; Interleukin 6; IL-6; IFNB2; BSF2; HGF; HSF; B-Cell Stimulatory Factor 2; CTL Differentiation Factor; Hybridoma Growth Factor; Interferon Beta-2; Interleukin-6; IFN-Beta-2; BSF-2; CDF; Interleukin 6 (Interferon, Beta 2); B-Cell Differentiation Factor; Interferon, Beta 2; Interleukin BSF-2

Background

Gene Name: IL6

NCBI Gene Entry: [3569](#)

UniProt Entry: [P05231](#)

Application Information

Clonality: Rat monoclonal antibody

Clone ID: 24GB15505

Species Reactivity: Human

Applications Tested: Intracellular Staining for Flow Cytometry ☐ ICFC ☐

Conjugation: PE

Immunogen

The immunogen information is proprietary.

Isotype

Rat IgG2a ☐ k

Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

Storage

Store at 4 °C and protected from prolonged exposure to light.

Recommended Dilutions

5 µL/Tests

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

