

Anti-CD8 Mouse Monoclonal Antibody, PE-Cy5.5 Conjugated



Catalog #: F0628-PC5.5

Aliases

CD8A; CD8 Subunit Alpha; T-Cell Surface Glycoprotein CD8 Alpha Chain; CD8alpha; P32; T-Lymphocyte Differentiation Antigen T8/Leu-2; CD8 Antigen, Alpha Polypeptide (P32); CD8a Molecule; CD8; Leu2 T-Lymphocyte Antigen; OKT8 T-Cell Antigen; T-Cell Antigen Leu2; T Cell Co-Receptor; T8 T-Cell Antigen; CD8a Antigen; IMD116; Leu2; MAL

Background

Gene Name: CD8A

NCBI Gene Entry: [925](#)

UniProt Entry: [P01732](#)

Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB15025

Species Reactivity: Human

Applications Tested: Flow Cytometry (FCM)

Conjugation: PE-Cy5.5

Immunogen

The immunogen information is proprietary.

Isotype

Mouse IgG1

Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

Storage

Store at 4 °C and protected from prolonged exposure to light.

Recommended Dilutions

5 µL/Tests

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041


ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-CD8 Mouse Monoclonal Antibody, PE-Cy5.5 Conjugated

PAGE 2



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM