

# Anti-CD11b Mouse Monoclonal Antibody, PerCP Conjugated



## Catalog #: F0641-PerCP

### Aliases

ITGAM; Integrin Subunit Alpha M; MAC-1; CD11b; CR3A; Integrin, Alpha M (Complement Component 3 Receptor 3 Subunit); Cell Surface Glycoprotein MAC-1 Subunit Alpha; Complement Component 3 Receptor 3 Subunit; Macrophage-1 Antigen Alpha Subunit; CD11 Antigen-Like Family Member B; Leukocyte Adhesion Receptor MO1; Integrin Alpha-M; CR-3 Alpha Chain; Integrin, Alpha M (Complement Component Receptor 3, Alpha; Also Known As CD11b (P170), Macrophage Antigen Alpha Polypeptide); Neutrophil Adherence Receptor Alpha-M Subunit; Macrophage Antigen Alpha Polypeptide; Neutrophil Adherence Receptor; Antigen CD11b (P170); CD11b Antigen; MAC1A; SLEB6; CD11B; MO1A

### Background

Gene Name: ITGAM

NCBI Gene Entry: [3684](#)

UniProt Entry: [P11215](#)

### Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB14465

Species Reactivity: Human

Applications Tested: Flow Cytometry (FCM)

Conjugation: PerCP

### Immunogen

The immunogen information is proprietary.

### Isotype

Mouse IgG1

### Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

### Storage

Store at 4 °C and protected from prolonged exposure to light.

### Recommended Dilutions

5 µL/Tests

---

#### SUPPORT

SUPPORT@GENUINBIOTECH.COM  
TEL: +1-540-855-7041

#### ORDERS

SALES@GENUINBIOTECH.COM  
FAX: +1-540-855-7041


[WWW.GENUINBIOTECH.COM](http://WWW.GENUINBIOTECH.COM)

# Anti-CD11b Mouse Monoclonal Antibody, PerCP Conjugated

PAGE 2

**Note:** This product is for research use only.

## Validation Data



### SUPPORT

SUPPORT@GENUINBIOTECH.COM  
TEL: +1-540-855-7041

### ORDERS

SALES@GENUINBIOTECH.COM  
FAX: +1-540-855-7041

[WWW.GENUINBIOTECH.COM](http://WWW.GENUINBIOTECH.COM)