

# Anti-CD61 Mouse Monoclonal Antibody, PerCP Conjugated



## Catalog #: F0661-PerCP

### Aliases

ITGB3; Integrin Subunit Beta 3; GPIIIa; CD61; GP3A; Integrin, Beta 3 (Platelet Glycoprotein IIIa, Antigen CD61); Platelet Membrane Glycoprotein IIIa; Integrin Beta-3; Antigen CD61; Integrin Beta Chain, Beta 3; Platelet Glycoprotein IIIa; Integrin Beta 3; CD61 Antigen; BDPLT16; BDPLT24; BDPLT2; GT2; GT

### Background

Gene Name: ITGB3

NCBI Gene Entry: [3690](#)

UniProt Entry: [P05106](#)

### Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB14720

Species Reactivity: Human

Applications Tested: Flow Cytometry (FCM)

Conjugation: PerCP

### Immunogen

The immunogen information is proprietary.

### Isotype

Mouse IgG1

### Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

### Storage

Store at 4 °C and protected from prolonged exposure to light.

### Recommended Dilutions

5 µL/Tests

**Note:** This product is for research use only.

### Validation Data

---

#### SUPPORT

SUPPORT@GENUINBIOTECH.COM  
TEL: +1-540-855-7041


#### ORDERS

SALES@GENUINBIOTECH.COM  
FAX: +1-540-855-7041

[WWW.GENUINBIOTECH.COM](http://WWW.GENUINBIOTECH.COM)

## Anti-CD61 Mouse Monoclonal Antibody, PerCP Conjugated

PAGE 2



### SUPPORT

SUPPORT@GENUINBIOTECH.COM  
TEL: +1-540-855-7041

### ORDERS

SALES@GENUINBIOTECH.COM  
FAX: +1-540-855-7041

[WWW.GENUINBIOTECH.COM](http://WWW.GENUINBIOTECH.COM)