

Anti-CD279 Mouse Monoclonal Antibody, PerCP Conjugated



Catalog #: F0662-PerCP

Aliases

PDCD1; Programmed Cell Death 1; PD1; CD279; HSLE1; PD-1; Systemic Lupus Erythematosus Susceptibility 2; Programmed Cell Death Protein 1; Protein PD-1; SLEB2; HPD-1; Programmed Cell Death 1 Protein; CD279 Antigen; AIMTBS; HPD-L

Background

Gene Name: PDCD1

NCBI Gene Entry: [5133](#)

UniProt Entry: [Q15116](#)

Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB14940

Species Reactivity: Human, african green, baboon, chimpanzee, marmoset, cynomolgus, rhesus, Saimiri sciureus

Applications Tested: Flow Cytometry (FCM)

Conjugation: PerCP

Immunogen

The immunogen information is proprietary.

Isotype

Mouse IgG1, κ

Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

Storage

Store at 4 °C and protected from prolonged exposure to light.

Recommended Dilutions

5 μ L/Tests

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041


ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-CD279 Mouse Monoclonal Antibody, PerCP Conjugated

PAGE 2



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM