

Anti-CD3 Mouse Monoclonal Antibody, PerCP-Cy5.5 Conjugated



Catalog #: F0665-PC5.5

Aliases

CD3G; CD3 Gamma Subunit Of T-Cell Receptor Complex; T-Cell Surface Glycoprotein CD3 Gamma Chain; CD3-GAMMA; CD3GAMMA; CD3g Antigen, Gamma Polypeptide (TiT3 Complex); CD3g Molecule, Gamma (CD3-TCR Complex); T-Cell Receptor T3 Gamma Chain; T3G; T-Cell Antigen Receptor Complex, Gamma Subunit Of T3; CD3g Molecule, Epsilon (CD3-TCR Complex); CD3g Molecule; CD3g Antigen; IMD17

Background

Gene Name: CD3G

NCBI Gene Entry: [917](#)

UniProt Entry: [P09693](#)

Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB14405

Species Reactivity: Human

Applications Tested: Flow Cytometry (FCM)

Conjugation: PerCP-Cy5.5

Immunogen

The immunogen information is proprietary.

Isotype

Mouse IgG2a

Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

Storage

Store at 4 °C and protected from prolonged exposure to light.

Recommended Dilutions

5 µL/Tests

Note: This product is for research use only.

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS


SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-CD3 Mouse Monoclonal Antibody, PerCP-Cy5.5 Conjugated

PAGE 2

Validation Data



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM