

Anti-CD4 Mouse Monoclonal Antibody, PerCP-Cy5.5 Conjugated



Catalog #: F0667-PC5.5

Aliases

CD4 □ CD4 Molecule □ T-Cell Surface Glycoprotein CD4 □ Leu-3 □ T4 □ T-Cell Surface Antigen T4/Leu-3 □ CD4 Antigen (P55) □ CD4 Receptor □ CD4 Antigen □ CD4mut □ IMD79 □ OKT4D

Background

Gene Name: CD4

NCBI Gene Entry: [920](#)

UniProt Entry: [P01730](#)

Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB15445

Species Reactivity: Human

Applications Tested: Flow Cytometry (FCM)

Conjugation: PerCP-Cy5.5

Immunogen

The immunogen information is proprietary.

Isotype

Mouse IgG1, κ

Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

Storage

Store at 4 °C and protected from prolonged exposure to light.

Recommended Dilutions

5 µL/Tests

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041


ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-CD4 Mouse Monoclonal Antibody, PerCP-Cy5.5 Conjugated

PAGE 2



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM