

Anti-CD5 Mouse Monoclonal Antibody, PerCP-Cy5.5 Conjugated



Catalog #: F0668-PC5.5

Aliases

CD5; CD5 Molecule; LEU1; T1; T-Cell Surface Glycoprotein CD5; Lymphocyte Antigen T1/Leu-1; CD5 Antigen (P56-62); Epididymis Secretory Sperm Binding Protein; CD5 Antigen

Background

Gene Name: CD5

NCBI Gene Entry: [921](#)

UniProt Entry: [P06127](#)

Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB14425

Species Reactivity: Human

Applications Tested: Flow Cytometry (FCM)

Conjugation: PerCP-Cy5.5

Immunogen

The immunogen information is proprietary.

Isotype

Mouse IgG1, κ

Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

Storage

Store at 4 °C and protected from prolonged exposure to light.

Recommended Dilutions

5 μ L/Tests

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041


ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-CD5 Mouse Monoclonal Antibody, PerCP-Cy5.5 Conjugated

PAGE 2



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM