

Anti-CD10 Mouse Monoclonal Antibody, PerCP-Cy5.5 Conjugated



Catalog #: F0670-PC5.5

Aliases

MME; Membrane Metalloendopeptidase; CALLA, NEP; Neprilysin; CD10; Common Acute Lymphocytic Leukemia Antigen; Neutral Endopeptidase 24.11; Skin Fibroblast Elastase; Neutral Endopeptidase; Atriopeptidase; Enkephalinase; EC 3.4.24.11; SFE; Membrane Metallo-Endopeptidase (Neutral Endopeptidase, Enkephalinase, CALLA, CD10); Neprilysin-390; Neprilysin-411; CD10 Antigen; EC 3.4.24 47; CMT2T; SCA43; EPN

Background

Gene Name: MME

NCBI Gene Entry: [4311](#)

UniProt Entry: [P08473](#)

Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB15535

Species Reactivity: Human

Applications Tested: Flow Cytometry (FCM)

Conjugation: PerCP-Cy5.5

Immunogen

The immunogen information is proprietary.

Isotype

Mouse IgG1, κ

Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

Storage

Store at 4 °C and protected from prolonged exposure to light.

Recommended Dilutions

5 μ L/Tests

Note: This product is for research use only.

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS


SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-CD10 Mouse Monoclonal Antibody, PerCP-Cy5.5 Conjugated

PAGE 2

Validation Data



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM