

Anti-D14 Mouse Monoclonal Antibody, PerCP-Cy5.5 Conjugated



Catalog #: F0672-PC5.5

Aliases

CD14; CD14 Molecule; Myeloid Cell-Specific Leucine-Rich Glycoprotein; Monocyte Differentiation Antigen CD14; CD14 Antigen; My23 Antigen

Background

Gene Name: CD14

NCBI Gene Entry: [929](#)

UniProt Entry: [P08571](#)

Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB15540

Species Reactivity: Human

Applications Tested: Flow Cytometry (FCM)

Conjugation: PerCP-Cy5.5

Immunogen

The immunogen information is proprietary.

Isotype

Mouse IgG1, κ

Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

Storage

Store at 4 °C and protected from prolonged exposure to light.

Recommended Dilutions

5 μ L/Tests

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041


ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-D14 Mouse Monoclonal Antibody, PerCP-Cy5.5 Conjugated

PAGE 2



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM