

# Anti-CD15 Mouse Monoclonal Antibody, PerCP-Cy5.5 Conjugated



**Catalog #: F0673-PC5.5**

## Aliases

FUT4; Fucosyltransferase 4; Galactoside 3-L-Fucosyltransferase; FUC-TIV; FCT3A; ELFT; 4-Galactosyl-N-Acetylglucosaminide 3-Alpha-L-Fucosyltransferase; Alpha (1,3) Fucosyltransferase, Myeloid-Specific; Alpha-(1,3)-Fucosyltransferase 4; ELAM-1 Ligand Fucosyltransferase; ELAM Ligand Fucosyltransferase; Fucosyltransferase IV; FucT-IV; CD15; Stage-Specific Embryonic Antigen 1; EC 2.4.1.152; EC 2.4.1.65; EC 2.4.1; Lewis X; Fuc-TIV; SSEA-1; FUTIV; LeX

## Background

Gene Name: FUT4

NCBI Gene Entry: [2526](#)

UniProt Entry: [P22083](#)

## Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB15545

Species Reactivity: Human, chimpanzee

Applications Tested: Flow Cytometry (FCM)

Conjugation: PerCP-Cy5.5

## Immunogen

The immunogen information is proprietary.

## Isotype

Mouse IgG1,  $\kappa$

## Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

## Storage

Store at 4 °C and protected from prolonged exposure to light.

## Recommended Dilutions

5  $\mu$ L/Tests

**Note:** This product is for research use only.

---

### SUPPORT

SUPPORT@GENUINBIOTECH.COM  
TEL: +1-540-855-7041

### ORDERS


SALES@GENUINBIOTECH.COM  
FAX: +1-540-855-7041

[WWW.GENUINBIOTECH.COM](http://WWW.GENUINBIOTECH.COM)

# Anti-CD15 Mouse Monoclonal Antibody, PerCP-Cy5.5 Conjugated

PAGE 2

## Validation Data



### SUPPORT

SUPPORT@GENUINBIOTECH.COM  
TEL: +1-540-855-7041

### ORDERS

SALES@GENUINBIOTECH.COM  
FAX: +1-540-855-7041

[WWW.GENUINBIOTECH.COM](http://WWW.GENUINBIOTECH.COM)