

Anti-CD20 Mouse Monoclonal Antibody, PerCP-Cy5.5 Conjugated



Catalog #: F0675-PC5.5

Aliases

MS4A1; Membrane Spanning 4-Domains A1; Bp35; FMC7; CD20; B1; Membrane-Spanning 4-Domains, Subfamily A, Member 1; Leukocyte Surface Antigen Leu-16; B-Lymphocyte Antigen CD20; CD20 Antigen; Membrane-Spanning 4-Domains Subfamily A Member 1; B-Lymphocyte Cell-Surface Antigen B1; B-Lymphocyte Surface Antigen B1; CD20 Receptor; LEU-16; CVID5; S7

Background

Gene Name: MS4A1

NCBI Gene Entry: [931](#)

UniProt Entry: [P11836](#)

Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB14525

Species Reactivity: Human

Applications Tested: Flow Cytometry (FCM)

Conjugation: PerCP-Cy5.5

Immunogen

The immunogen information is proprietary.

Isotype

Mouse IgG2b, κ

Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

Storage

Store at 4 °C and protected from prolonged exposure to light.

Recommended Dilutions

5 μ L/Tests

Note: This product is for research use only.

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS


SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-CD20 Mouse Monoclonal Antibody, PerCP-Cy5.5 Conjugated

PAGE 2

Validation Data



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM