

# Anti-CD23 Mouse Monoclonal Antibody, PerCP-Cy5.5 Conjugated



**Catalog #: F0676-PC5.5**

## Aliases

FCER2; Fc Epsilon Receptor II; CLEC4J; FcepsilonRII; FCErII; CD23A; CD23; FCE2; Fc Fragment Of IgE, Low Affinity II, Receptor For (CD23); Low Affinity Immunoglobulin Epsilon Fc Receptor; C-Type Lectin Domain Family 4 Member J; Immunoglobulin E-Binding Factor; Fc Fragment Of IgE Receptor II; Lymphocyte IgE Receptor; Fc-Epsilon-RII; CD23 Antigen; BLAST-2; IGEBF; Fc Fragment Of IgE, Low Affinity II, Receptor For (CD23A); C-Type Lectin Domain Family 4, Member J; Immunoglobulin Epsilon-Chain; Fc Epsilon RII

## Background

Gene Name: FCER2

NCBI Gene Entry: [2208](#)

UniProt Entry: [P06734](#)

## Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB14540

Species Reactivity: Human

Applications Tested: Flow Cytometry (FCM)

Conjugation: PerCP-Cy5.5

## Immunogen

The immunogen information is proprietary.

## Isotype

Mouse IgG1,  $\kappa$

## Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

## Storage

Store at 4 °C and protected from prolonged exposure to light.

## Recommended Dilutions

5  $\mu$ L/Tests

---

### SUPPORT

SUPPORT@GENUINBIOTECH.COM  
TEL: +1-540-855-7041

### ORDERS

SALES@GENUINBIOTECH.COM  
FAX: +1-540-855-7041


[WWW.GENUINBIOTECH.COM](http://WWW.GENUINBIOTECH.COM)

## Anti-CD23 Mouse Monoclonal Antibody, PerCP-Cy5.5 Conjugated

PAGE 2

**Note:** This product is for research use only.

### Validation Data



#### SUPPORT

SUPPORT@GENUINBIOTECH.COM  
TEL: +1-540-855-7041

#### ORDERS

SALES@GENUINBIOTECH.COM  
FAX: +1-540-855-7041

[WWW.GENUINBIOTECH.COM](http://WWW.GENUINBIOTECH.COM)