

Anti-CD28 Mouse Monoclonal Antibody, PerCP-Cy5.5 Conjugated



Catalog #: F0679-PC5.5

Aliases

CD28; CD28 Molecule; T-Cell-Specific Surface Glycoprotein CD28; T-Cell-Specific Surface Glycoprotein; CD28 Antigen (Tp44); CD28 Antigen; IMD123; Tp44; TP44

Background

Gene Name: CD28

NCBI Gene Entry: [940](#)

UniProt Entry: [P10747](#)

Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB14575

Species Reactivity: Human, baboon, capuchin monkey, chimpanzee, cynomolgus, rhesus, saimiri sciureus

Applications Tested: Flow Cytometry (FCM)

Conjugation: PerCP-Cy5.5

Immunogen

The immunogen information is proprietary.

Isotype

Mouse IgG1, κ

Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

Storage

Store at 4 °C and protected from prolonged exposure to light.

Recommended Dilutions

5 μ L/Tests

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041


ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-CD28 Mouse Monoclonal Antibody, PerCP-Cy5.5 Conjugated

PAGE 2



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM