

Anti-CD33 Mouse Monoclonal Antibody, PerCP-Cy5.5 Conjugated



Catalog #: F0680-PC5.5

Aliases

CD33; CD33 Molecule; SIGLEC3; CD33rSiglec; SIGLEC-3; P67; ialic Acid-Binding Ig-Like Lectin 3; Myeloid Cell Surface Antigen CD33; CD33 Antigen (Gp67); FLJ00391; Gp67; Sialic Acid Binding Ig-Like Lectin 3; CD33 Molecule Transcript; CD33 Antigen; Siglec-3

Background

Gene Name: CD33

NCBI Gene Entry: [945](#)

UniProt Entry: [P20138](#)

Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB15390

Species Reactivity: Human

Applications Tested: Flow Cytometry (FCM)

Conjugation: PerCP-Cy5.5

Immunogen

The immunogen information is proprietary.

Isotype

Mouse IgG1, κ

Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

Storage

Store at 4 °C and protected from prolonged exposure to light.

Recommended Dilutions

5 μ L/Tests

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041


ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-CD33 Mouse Monoclonal Antibody, PerCP-Cy5.5 Conjugated

PAGE 2



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM