

Anti-CD38 Mouse Monoclonal Antibody, PerCP-Cy5.5 Conjugated



Catalog #: F0682-PC5.5

Aliases

CD38; CD38 Molecule; ADP-Ribosyl Cyclase 1; CADPR1; ADP-Ribosyl Cyclase/Cyclic ADP-Ribose Hydrolase 1; 2'-Phospho-Cyclic-ADP-Ribose Transferase; 2'-Phospho-ADP-Ribosyl Cyclase; Cyclic ADP-Ribose Hydrolase 1; NAD(+) Nucleosidase; CD38 Antigen (P45); ADPRC 1; 2'-Phospho-ADP-Ribosyl Cyclase/2'-Phospho-Cyclic-ADP-Ribose Transferase; Ecto-Nicotinamide Adenine Dinucleotide Glycohydrolase; Cluster Of Differentiation 38; CADPR Hydrolase 1; CD38 Antigen; EC 2.4.99.20; EC 3.2.2.-; EC 3.2.2.6; ADPRC1; T10

Background

Gene Name: CD38

NCBI Gene Entry: [952](#)

UniProt Entry: [P28907](#)

Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB14615

Species Reactivity: Human, chimpanzee, horse, bovine

Applications Tested: Flow Cytometry (FCM)

Conjugation: PerCP-Cy5.5

Immunogen

The immunogen information is proprietary.

Isotype

Mouse IgG1

Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

Storage

Store at 4 °C and protected from prolonged exposure to light.

Recommended Dilutions

5 µL/Tests

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041


WWW.GENUINBIOTECH.COM

Anti-CD38 Mouse Monoclonal Antibody, PerCP-Cy5.5 Conjugated

PAGE 2

Note: This product is for research use only.

Validation Data



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM