

Anti-CD56 Mouse Monoclonal Antibody, PerCP-Cy5.5 Conjugated



Catalog #: F0685-PC5.5

Aliases

NCAM1; Neural Cell Adhesion Molecule 1; NCAM; CD56; Antigen Recognized By Monoclonal Antibody 5.1H11; Neural Cell Adhesion Molecule, NCAM; CD56 Antigen; N-CAM-1; NCAM-1; MSK39

Background

Gene Name: NCAM1

NCBI Gene Entry: [4684](#)

UniProt Entry: [P13591](#)

Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB15555

Species Reactivity: Human

Applications Tested: Flow Cytometry (FCM)

Conjugation: PerCP-Cy5.5

Immunogen

The immunogen information is proprietary.

Isotype

Mouse IgG2a

Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

Storage

Store at 4 °C and protected from prolonged exposure to light.

Recommended Dilutions

5 µL/Tests

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041


ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-CD56 Mouse Monoclonal Antibody, PerCP-Cy5.5 Conjugated

PAGE 2



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM