

Anti-CD64 Mouse Monoclonal Antibody, PerCP-Cy5.5 Conjugated



Catalog #: F0686-PC5.5

Aliases

FCGR1A; Fc Gamma Receptor Ia; FcgammaRIa; FCGR1; FCG1; FcgammaRI; CD64A; CD64; Fc Fragment Of IgG, High Affinity Ia, Receptor For (CD64); Fc Fragment Of IgG, High Affinity Ia, Receptor (CD64); High Affinity Immunoglobulin Gamma Fc Receptor I; Fc Fragment Of IgG Receptor Ia; IgG Fc Receptor I; Fc-Gamma RIA; Fc-Gamma RI; IGFR1; Fc-Gamma Receptor I A1; CD64 Antigen; Fc Gamma RI; FCRI; FcRI

Background

Gene Name: FCGR1A

NCBI Gene Entry: [2209](#)

UniProt Entry: [P12314](#)

Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB14735

Species Reactivity: Human, baboon, capuchin monkey, chimpanzee, cynomolgus, rhesus, saimiri sciureus

Applications Tested: Flow Cytometry (FCM)

Conjugation: PerCP-Cy5.5

Immunogen

The immunogen information is proprietary.

Isotype

Mouse IgG1

Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

Storage

Store at 4 °C and protected from prolonged exposure to light.

Recommended Dilutions

5 µL/Tests

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041


WWW.GENUINBIOTECH.COM

Anti-CD64 Mouse Monoclonal Antibody, PerCP-Cy5.5 Conjugated

PAGE 2

Note: This product is for research use only.

Validation Data



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM