

Anti-CD117 Mouse Monoclonal Antibody, PerCP-Cy5.5 Conjugated



Catalog #: F0688-PC5.5

Aliases

KIT; KIT Proto-Oncogene, Receptor Tyrosine Kinase; SCFR; V-Kit Hardy-Zuckerman 4 Feline Sarcoma Viral Oncogene Homolog; Mast/Stem Cell Growth Factor Receptor Kit; CD117; C-Kit; PBT; Tyrosine-Protein Kinase Kit; Piebald Trait Protein; Proto-Oncogene C-Kit; EC 2.7.10.1; P145 C-Kit; V-Kit Hardy-Zuckerman 4 Feline Sarcoma Viral Oncogene-Like Protein; Proto-Oncogene Tyrosine-Protein Kinase Kit; C-Kit Protooncogene; Piebald Trait; CD117 Antigen; EC 2.7.10; MASTC

Background

Gene Name: KIT

NCBI Gene Entry: [3815](#)

UniProt Entry: [P10721](#)

Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB14830

Species Reactivity: Human

Applications Tested: Flow Cytometry (FCM)

Conjugation: PerCP-Cy5.5

Immunogen

The immunogen information is proprietary.

Isotype

Mouse IgG1, κ

Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

Storage

Store at 4 °C and protected from prolonged exposure to light.

Recommended Dilutions

5 μ L/Tests

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041


WWW.GENUINBIOTECH.COM

Anti-CD117 Mouse Monoclonal Antibody, PerCP-Cy5.5 Conjugated

PAGE 2

Note: This product is for research use only.

Validation Data



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM