

Anti-CD123 Mouse Monoclonal Antibody, PerCP-Cy5.5 Conjugated



Catalog #: F0689-PC5.5

Aliases

IL3RA; Interleukin 3 Receptor Subunit Alpha; CD123; Interleukin 3 Receptor, Alpha (Low Affinity); Interleukin-3 Receptor Subunit Alpha; IL-3 Receptor Subunit Alpha; IL-3R Subunit Alpha; CD123 Antigen; IL-3R-Alpha; IL3R; HIL-3Ra; IL3RAY; IL-3RA; IL3RX; IL3RY

Background

Gene Name: IL3RA

NCBI Gene Entry: [3563](#)

UniProt Entry: [P26951](#)

Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB14835

Species Reactivity: Human

Applications Tested: Flow Cytometry (FCM)

Conjugation: PerCP-Cy5.5

Immunogen

The immunogen information is proprietary.

Isotype

Mouse IgG1, κ

Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

Storage

Store at 4 °C and protected from prolonged exposure to light.

Recommended Dilutions

5 μ L/Tests

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041


ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-CD123 Mouse Monoclonal Antibody, PerCP-Cy5.5 Conjugated

PAGE 2



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM