

Anti-CD127 Mouse Monoclonal Antibody, PerCP-Cy5.5 Conjugated



Catalog #: F0690-PC5.5

Aliases

IL7R; Interleukin 7 Receptor; CDw127; Interleukin-7 Receptor Subunit Alpha; IL-7Ralpha; IL7Ralpha; Lnc-IL7R; SIL-7R; CD127; IL7RA; Soluble Interleukin-7 Receptor; IL-7 Receptor Subunit Alpha; IL-7R Subunit Alpha; CD127 Antigen; IL-7R-Alpha; Interleukin 7 Receptor Alpha Chain; IMD104; IL-7RA; ILRA

Background

Gene Name: IL7R

NCBI Gene Entry: [3575](#)

UniProt Entry: [P16871](#)

Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB14840

Species Reactivity: Human

Applications Tested: Flow Cytometry (FCM)

Conjugation: PerCP-Cy5.5

Immunogen

The immunogen information is proprietary.

Isotype

Mouse IgG1, κ

Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

Storage

Store at 4 °C and protected from prolonged exposure to light.

Recommended Dilutions

5 μ L/Tests

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041


ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-CD127 Mouse Monoclonal Antibody, PerCP-Cy5.5 Conjugated

PAGE 2



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM