

Anti-CD9 Mouse Monoclonal Antibody, Pacific Blue Conjugated



Catalog #: F0698-PB

Aliases

CD9; CD9 Molecule; TSPAN29; MRP-1; Tetraspanin-29; MIC3; P24; Cell Growth-Inhibiting Gene 2 Protein; Motility Related Protein-1; Leukocyte Antigen MIC3; CD9 Antigen (P24); CD9 Antigen; 5H9 Antigen; BA2; Motility-Related Protein; BA-2/P24 Antigen; Antigen CD9; TSPAN-29; Tspan-29; DRAP-27; BTCC-1

Background

Gene Name: CD9

NCBI Gene Entry: [928](#)

UniProt Entry: [P21926](#)

Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB14450

Species Reactivity: Human

Applications Tested: Flow Cytometry (FCM)

Conjugation: Pacific Blue

Immunogen

The immunogen information is proprietary.

Isotype

Mouse IgG1, κ

Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

Storage

Store at 4 °C and protected from prolonged exposure to light.

Recommended Dilutions

5 μ L/Tests

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041


ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-CD9 Mouse Monoclonal Antibody, Pacific Blue Conjugated

PAGE 2



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM