

Anti-CD18 Mouse Monoclonal Antibody, Pacific Blue Conjugated



Catalog #: F0711-PB

Aliases

ITGB2; Integrin Subunit Beta 2; LFA-1; MAC-1; CD18; MFI7; Integrin, Beta 2 (Complement Component 3 Receptor 3 And 4 Subunit); Complement Component 3 Receptor 3 And 4 Subunit; Integrin Beta-2; Integrin, Beta 2 (Antigen CD18 (P95), Lymphocyte Function-Associated Antigen 1; Macrophage Antigen 1 (Mac-1) Beta Subunit); Cell Surface Adhesion Glycoprotein (LFA-1/CR3/P150,959 Beta Subunit Precursor); Cell Surface Adhesion Glycoproteins LFA-1/CR3/P150,95 Subunit Beta; Leukocyte-Associated Antigens CD18/11A, CD18/11B, CD18/11C; Leukocyte Cell Adhesion Molecule CD18; Complement Receptor C3 Beta-Subunit; Complement Receptor C3 Subunit Beta; Integrin Beta Chain, Beta 2; D18 Antigen; LCAMB; MF17; LAD

Background

Gene Name: ITGB2

NCBI Gene Entry: [3689](#)

UniProt Entry: [P05107](#)

Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB15565

Species Reactivity: Human, pig

Applications Tested: Flow Cytometry (FCM)

Conjugation: Pacific Blue

Immunogen

The immunogen information is proprietary.

Isotype

Mouse IgG1 κ

Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

Storage

Store at 4 °C and protected from prolonged exposure to light.

Recommended Dilutions

5 μ L/Tests

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041


WWW.GENUINBIOTECH.COM

Anti-CD18 Mouse Monoclonal Antibody, Pacific Blue Conjugated

PAGE 2

Note: This product is for research use only.

Validation Data



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM