

Anti-CD21 Mouse Monoclonal Antibody, Pacific Blue Conjugated



Catalog #: F0714-PB

Aliases

CR2; Complement C3d Receptor 2; C3DR; CD21; Complement Component (3d/Epstein Barr Virus) Receptor 2; Complement Component 3d Receptor 2; Epstein-Barr Virus Receptor; Complement Receptor Type 2; EBV Receptor; Complement C3d Receptor; CD21 Antigen; CVID7; SLEB9; Cr2; CR

Background

Gene Name: CR2

NCBI Gene Entry: [1380](#)

UniProt Entry: [P20023](#)

Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB14530

Species Reactivity: Human

Applications Tested: Flow Cytometry (FCM)

Conjugation: Pacific Blue

Immunogen

The immunogen information is proprietary.

Isotype

Mouse IgG1

Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

Storage

Store at 4 °C and protected from prolonged exposure to light.

Recommended Dilutions

5 µL/Tests

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041


ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-CD21 Mouse Monoclonal Antibody, Pacific Blue Conjugated

PAGE 2



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM