

Anti-CD26 Mouse Monoclonal Antibody, Pacific Blue Conjugated



Catalog #: F0718-PB

Aliases

DPP4; Dipeptidyl Peptidase 4; Adenosine Deaminase Complexing Protein 2; DPPIV; ADCP2; CD26; Dipeptidylpeptidase IV (CD26, Adenosine Deaminase Complexing Protein 2); T-Cell Activation Antigen CD26; Dipeptidyl Peptidase IV; EC 3.4.14.5; ADCP-2; DPP IV; ADABP; TP103; Post-Proline Dipeptidyl Aminopeptidase IV; Xaa-Pro-Dipeptidylaminopeptidase; Gly-Pro Naphthylamidase; Dipeptidyl-Peptidase 4; Dipeptidylpeptidase 4; CD26 Antigen

Background

Gene Name: DPP4

NCBI Gene Entry: [1803](#)

UniProt Entry: [P27487](#)

Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB15130

Species Reactivity: Human

Applications Tested: Flow Cytometry (FCM)

Conjugation: Pacific Blue

Immunogen

The immunogen information is proprietary.

Isotype

Mouse IgG2a, κ

Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

Storage

Store at 4 °C and protected from prolonged exposure to light.

Recommended Dilutions

5 μ L/Tests

Note: This product is for research use only.

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS


SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-CD26 Mouse Monoclonal Antibody, Pacific Blue Conjugated

PAGE 2

Validation Data



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM