

Anti-CD32 Mouse Monoclonal Antibody, Pacific Blue Conjugated



Catalog #: F0721-PB

Aliases

FCGR2A; Fc Gamma Receptor IIa ; Fc-Gamma-RIIa; IGFR2; CDw32; CD32; FcgammaRIIa; FCGR2A1; CD32A; FCG2; Low Affinity Immunoglobulin Gamma Fc Region Receptor II-A; Fc Fragment Of IgG, Low Affinity IIa, Receptor (CD32); Immunoglobulin G Fc Receptor II; Fc Fragment Of IgG Receptor IIa; IgG Fc Receptor II-A; FcRII-A; FCGR2; Fc Fragment Of IgG, Low Affinity IIa, Receptor For (CD32); Fc Gamma Receptor RIIa3; Fc-Gamma RII-A; CD32 Antigen; FcGR

Background

Gene Name: FCGR2A

NCBI Gene Entry: [2212](#)

UniProt Entry: [P12318](#)

Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB14595

Species Reactivity: Human

Applications Tested: Flow Cytometry (FCM)

Conjugation: Pacific Blue

Immunogen

The immunogen information is proprietary.

Isotype

Mouse IgG2b κ

Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

Storage

Store at 4 °C and protected from prolonged exposure to light.

Recommended Dilutions

5 μ L/Tests

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041


WWW.GENUINBIOTECH.COM

Anti-CD32 Mouse Monoclonal Antibody, Pacific Blue Conjugated

PAGE 2

Note: This product is for research use only.

Validation Data



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM