

Anti-CD44 Mouse Monoclonal Antibody, Pacific Blue Conjugated



Catalog #: F0728-PB

Aliases

CD44; CD44 Molecule (IN Blood Group); HUTCH-I; CDw44; Hermes-1; HUTCH-1; ECM-III; PGP-1; HCELL; CSPG8; H-CAM; MC56; Pgp1; MDU2; MDU3; MIC4; IN; Hematopoietic Cell E- And L-Selectin Ligand; GP90 Lymphocyte Homing/Adhesion Receptor; Chondroitin Sulfate Proteoglycan 8; Extracellular Matrix Receptor III; Homing Cell Adhesion Molecule; Heparan Sulfate Proteoglycan; Phagocytic Glycoprotein 1; Phagocyte Glycoprotein 1; Hyaluronate Receptor; In(Lu) Related-P80; Hermes Antigen; CD44 Antigen; ECMR-III; Epican; CD44R; LHR; CD44 Antigen (Homing Function And Indian Blood Group System); Homing Function And Indian Blood Group System; CD44 Molecule (Indian Blood Group); Extracellular Matrix Receptor-III; Cell Surface Glycoprotein CD44; Indian Blood Group Antigen; Phagocytic Glycoprotein I; Soluble CD44; PGP-I

Background

Gene Name: CD44

NCBI Gene Entry: [960](#)

UniProt Entry: [P16070](#)

Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB14645

Species Reactivity: Human

Applications Tested: Flow Cytometry (FCM)

Conjugation: Pacific Blue

Immunogen

The immunogen information is proprietary.

Isotype

Mouse IgG2a

Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

Storage

Store at 4 °C and protected from prolonged exposure to light.

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-CD44 Mouse Monoclonal Antibody, Pacific Blue Conjugated


PAGE 2

Recommended Dilutions

5 μ L/Tests

Note: This product is for research use only.

Validation Data



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM