

Anti-CD57 Mouse Monoclonal Antibody, Pacific Blue Conjugated



Catalog #: F0733-PB

Aliases

B3GAT1; Beta-1,3-Glucuronyltransferase 1; Galactosylgalactosylxylosylprotein 3-Beta-Glucuronosyltransferase 1; Glucuronosyltransferase P; GlcAT-P; UDP-GlcUA:Glycoprotein Beta-1,3-Glucuronyltransferase; EC 2.4.1.135; GLCATP; CD57; LEU7; CD57 Antigen; LEU7 Antigen; GlcUAT-P; GLCUATP; HNK1; NK-1; NK1

Background

Gene Name: B3GAT1

NCBI Gene Entry: [27087](#)

UniProt Entry: [Q9P2W7](#)

Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB14710

Species Reactivity: Human

Applications Tested: Flow Cytometry (FCM)

Conjugation: Pacific Blue

Immunogen

The immunogen information is proprietary.

Isotype

Mouse IgM

Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

Storage

Store at 4 °C and protected from prolonged exposure to light.

Recommended Dilutions

5 µL/Tests

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041


ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-CD57 Mouse Monoclonal Antibody, Pacific Blue Conjugated

PAGE 2



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM