

Anti-CD81 Mouse Monoclonal Antibody, Pacific Blue Conjugated



Catalog #: F0736-PB

Aliases

CD81; CD81 Molecule; TSPAN28; Tetraspanin-28; TAPA1; S5.7; CD81 Antigen (Target Of Antiproliferative Antibody 1); 26 KDa Cell Surface Protein TAPA-1; CD81 Antigen; Tspan-28; TAPA-1; Target Of The Antiproliferative Antibody 1; CVID6

Background

Gene Name: CD81

NCBI Gene Entry: [975](#)

UniProt Entry: [P60033](#)

Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB14770

Species Reactivity: Human

Applications Tested: Flow Cytometry (FCM)

Conjugation: Pacific Blue

Immunogen

The immunogen information is proprietary.

Isotype

Mouse IgG1

Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

Storage

Store at 4 °C and protected from prolonged exposure to light.

Recommended Dilutions

5 µL/Tests

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041


ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-CD81 Mouse Monoclonal Antibody, Pacific Blue Conjugated

PAGE 2



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM