

# Anti-CD278 Armenian Hamstere Monoclonal Antibody, Pacific Blue Conjugated



**Catalog #: F0740-PB**

## Aliases

ICOS; Inducible T Cell Costimulator; AILIM; Activation-Inducible Lymphocyte Immunomediatory Molecule; Inducible T-Cell Costimulator; CD278; Inducible T-Cell Co-Stimulator; Inducible Costimulator; CD278 Antigen; CVID1

## Background

Gene Name: ICOS

NCBI Gene Entry: [29851](#)

UniProt Entry: [Q9Y6W8](#)

## Application Information

Clonality: Armenian Hamstere monoclonal antibody

Clone ID: 24GB14935

Species Reactivity: Human

Applications Tested: Flow Cytometry (FCM)

Conjugation: Pacific Blue

## Immunogen

The immunogen information is proprietary.

## Isotype

Armenian Hamster IgG

## Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

## Storage

Store at 4 °C and protected from prolonged exposure to light.

## Recommended Dilutions

5 µL/Tests

**Note:** This product is for research use only.

## Validation Data

---

### SUPPORT

SUPPORT@GENUINBIOTECH.COM  
TEL: +1-540-855-7041


### ORDERS

SALES@GENUINBIOTECH.COM  
FAX: +1-540-855-7041

[WWW.GENUINBIOTECH.COM](http://WWW.GENUINBIOTECH.COM)

## Anti-CD278 Armenian Hamstere Monoclonal Antibody, Pacific Blue Conjugated

PAGE 2



### SUPPORT

SUPPORT@GENUINBIOTECH.COM  
TEL: +1-540-855-7041

### ORDERS

SALES@GENUINBIOTECH.COM  
FAX: +1-540-855-7041

[WWW.GENUINBIOTECH.COM](http://WWW.GENUINBIOTECH.COM)