

Anti-CD24 Rat Monoclonal Antibody, APC Conjugated



Catalog #: F0744-APC

Aliases

CD24; CD24 Molecule; CD24A; CD24 Antigen (Small Cell Lung Carcinoma Cluster 4 Antigen); Signal Transducer CD24; Small Cell Lung Carcinoma Cluster 4 Antigen; CD24 Antigen

Background

Gene Name: CD24

NCBI Gene Entry: [100133941](#)

UniProt Entry: [P25063](#)

Application Information

Clonality: Rat monoclonal antibody

Clone ID: 24GB15580

Species Reactivity: Mouse

Applications Tested: Flow Cytometry (FCM)

Conjugation: APC

Immunogen

The immunogen information is proprietary.

Isotype

Rat IgG2c, κ

Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

Storage

Store at 4 °C and protected from prolonged exposure to light.

Recommended Dilutions

5 μ L/Tests

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041


ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-CD24 Rat Monoclonal Antibody, APC Conjugated

PAGE 2



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM