

Anti-F4/80 Rat Monoclonal Antibody, APC Conjugated



Catalog #: F0749-APC

Aliases

ADGRE1; Adhesion G Protein-Coupled Receptor E1; TM7LN3; EMR1; Egf-Like Module Containing, Mucin-Like, Hormone Receptor-Like Sequence 1; Egf-Like Module Containing, Mucin-Like, Hormone Receptor-Like 1; EGF-Like Module-Containing Mucin-Like Hormone Receptor-Like 1; EGF-Like Module Receptor 1; EMR1 Hormone Receptor

Background

Gene Name: ADGRE1

NCBI Gene Entry: [2015](#)

UniProt Entry: [Q14246](#)

Application Information

Clonality: Rat monoclonal antibody

Clone ID: 24GB15605

Species Reactivity: Mouse

Applications Tested: Flow Cytometry (FCM)

Conjugation: APC

Immunogen

The immunogen information is proprietary.

Isotype

Rat IgG2a, κ

Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

Storage

Store at 4 °C and protected from prolonged exposure to light.

Recommended Dilutions

5 μ L/Tests

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041


ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-F4/80 Rat Monoclonal Antibody, APC Conjugated

PAGE 2



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM