

# Anti-I-A/I-E Rat Monoclonal Antibody, APC Conjugated



**Catalog #: F0752-APC**

## Aliases

IAb; Ia2; Ia-2; Abeta; H-2Ab; H2-Ab; Rmcs1; I-Abeta; Eb; Ia4; H2Eb; Ia-4; H-2Eb

## Background

Gene Name: H2-Ab1/H2-Eb1

NCBI Gene Entry: [14961/14969](#)

## Application Information

Clonality: Rat monoclonal antibody

Clone ID: 24GB15620

Species Reactivity: Mouse

Applications Tested: Flow Cytometry (FCM)

Conjugation: APC

## Immunogen

The immunogen information is proprietary.

## Isotype

Rat IgG2b,  $\kappa$

## Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

## Storage

Store at 4 °C and protected from prolonged exposure to light.

## Recommended Dilutions

5  $\mu$ L/Tests

**Note:** This product is for research use only.

## Validation Data

---

### SUPPORT

SUPPORT@GENUINBIOTECH.COM  
TEL: +1-540-855-7041


### ORDERS

SALES@GENUINBIOTECH.COM  
FAX: +1-540-855-7041

[WWW.GENUINBIOTECH.COM](http://WWW.GENUINBIOTECH.COM)

## Anti-I-A/I-E Rat Monoclonal Antibody, APC Conjugated

PAGE 2



### SUPPORT

SUPPORT@GENUINBIOTECH.COM  
TEL: +1-540-855-7041

### ORDERS

SALES@GENUINBIOTECH.COM  
FAX: +1-540-855-7041

[WWW.GENUINBIOTECH.COM](http://WWW.GENUINBIOTECH.COM)