

# Anti-LPAM-1 Rat Monoclonal Antibody, APC Conjugated



## Catalog #: F0753-APC

### Aliases

ITGA4; Integrin Subunit Alpha 4; CD49d; Alpha 4 Subunit Of VLA-4 Receptor; CD49 Antigen-Like Family Member D; VLA-4 Subunit Alpha; Integrin Alpha-IV; Integrin Alpha-4; Integrin, Alpha 4 (Antigen CD49D, Alpha 4 Subunit Of VLA-4 Receptor); Very Late Activation Protein 4 Receptor, Alpha 4 Subunit; Antigen CD49D, Alpha-4 Subunit Of VLA-4 Receptor; 269C Wild Type; Antigen CD49D; CD49d Antigen; CD49D; IA4; ITGB7; Integrin Subunit Beta 7; Gut Homing Receptor Beta Subunit; Integrin Beta-7; Integrin Beta 7 Subunit; Integrin, Beta 7

### Background

Gene Name: ITGA4/ITGB7

NCBI Gene Entry: [3676/3695](#)

UniProt Entry: [P13612/P26010](#)

### Application Information

Clonality: Rat monoclonal antibody

Clone ID: 24GB15625

Species Reactivity: Mouse

Applications Tested: Flow Cytometry (FCM)

Conjugation: APC

### Immunogen

The immunogen information is proprietary.

### Isotype

Rat IgG2a

### Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

### Storage

Store at 4 °C and protected from prolonged exposure to light.

### Recommended Dilutions

5 µL/Tests

---

#### SUPPORT

SUPPORT@GENUINBIOTECH.COM  
TEL: +1-540-855-7041

#### ORDERS

SALES@GENUINBIOTECH.COM  
FAX:+1-540-855-7041

WWW.GENUINBIOTECH.COM

# Anti-LPAM-1 Rat Monoclonal Antibody, APC Conjugated

**Note:** This product is for research use only.

## Validation Data

