

Anti-CD30 Armenian Hamster Monoclonal Antibody, FITC-Cy5 Conjugated



Catalog #: F0763-FITC

Aliases

TNFRSF8; TNF Receptor Superfamily Member 8; D1S166E; KI-1; CD30; Tumor Necrosis Factor Receptor Superfamily Member 8; Lymphocyte Activation Antigen CD30; CD30L Receptor; Ki-1 Antigen; Tumor Necrosis Factor Receptor Superfamily, Member 8; Cytokine Receptor CD30; CD30 Antigen

Background

Gene Name: TNFRSF8

NCBI Gene Entry: [943](#)

UniProt Entry: [P28908](#)

Application Information

Clonality: Armenian Hamster monoclonal antibody

Clone ID: 24GB15655

Species Reactivity: Mouse

Applications Tested: Flow Cytometry (FCM)

Conjugation: FITC-Cy5

Immunogen

The immunogen information is proprietary.

Isotype

Armenian Hamster IgG1

Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

Storage

Store at 4 °C and protected from prolonged exposure to light.

Recommended Dilutions

5 µL/Tests

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041


ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-CD30 Armenian Hamster Monoclonal Antibody, FITC-Cy5 Conjugated

PAGE 2



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM