

Anti-CD71 Rat Monoclonal Antibody, FITC-Cy5 Conjugated



Catalog #: F0765-FITC

Aliases

TFRC; Transferrin Receptor; P90; CD71; TFR1; Transferrin Receptor Protein 1; T9; TR; Transferrin Receptor (P90, CD71); CD71 Antigen; IMD46; TRFR; TfR1; Trfr; TFR; TfR

Background

Gene Name: TFRC

NCBI Gene Entry: [7037](#)

UniProt Entry: [P02786](#)

Application Information

Clonality: Rat monoclonal antibody

Clone ID: 24GB15660

Species Reactivity: Mouse

Applications Tested: Flow Cytometry (FCM)

Conjugation: FITC-Cy5

Immunogen

The immunogen information is proprietary.

Isotype

Rat IgG2a

Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

Storage

Store at 4 °C and protected from prolonged exposure to light.

Recommended Dilutions

5 µL/Tests

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041


ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-CD71 Rat Monoclonal Antibody, FITC-Cy5 Conjugated

PAGE 2



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM