

# Anti-CD154 Armenian Hamster Monoclonal Antibody, FITC-Cy5 Conjugated



## Catalog #: F0767-FITC

### Aliases

CD40LG; CD40 Ligand; CD40-L; CD40L; TRAP; TNF-Related Activation Protein; HCD40L; TNFSF5; CD154; HIGM1; T-BAM; Gp39; Tumor Necrosis Factor (Ligand) Superfamily Member 5; T-B Cell-Activating Molecule; CD40 Antigen Ligand; T-Cell Antigen Gp39; IMD3; Tumor Necrosis Factor (Ligand) Superfamily, Member 5 (Hyper-IgM Syndrome); Tumor Necrosis Factor Ligand Superfamily Member 5; Hyper-IgM Syndrome; CD154 Antigen; IGM

### Background

Gene Name: CD40LG

NCBI Gene Entry: [959](#)

UniProt Entry: [P29965](#)

### Application Information

Clonality: Armenian Hamster monoclonal antibody

Clone ID: 24GB15670

Species Reactivity: Mouse

Applications Tested: Flow Cytometry (FCM)

Conjugation: FITC-Cy5

### Immunogen

The immunogen information is proprietary.

### Isotype

Armenian Hamster IgG3

### Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

### Storage

Store at 4 °C and protected from prolonged exposure to light.

### Recommended Dilutions

5 µL/Tests

**Note:** This product is for research use only.

---

#### SUPPORT

SUPPORT@GENUINBIOTECH.COM  
TEL: +1-540-855-7041

#### ORDERS


SALES@GENUINBIOTECH.COM  
FAX: +1-540-855-7041

[WWW.GENUINBIOTECH.COM](http://WWW.GENUINBIOTECH.COM)

# Anti-CD154 Armenian Hamster Monoclonal Antibody, FITC-Cy5 Conjugated

PAGE 2

## Validation Data



### SUPPORT

SUPPORT@GENUINBIOTECH.COM  
TEL: +1-540-855-7041

### ORDERS

SALES@GENUINBIOTECH.COM  
FAX: +1-540-855-7041

[WWW.GENUINBIOTECH.COM](http://WWW.GENUINBIOTECH.COM)