

Anti-CD183 Armenian Hamster Monoclonal Antibody, FITC-Cy5 Conjugated



Catalog #: F0769-FITC

Aliases

CXCR3; C-X-C Motif Chemokine Receptor 3; CKR-L2; G Protein-Coupled Receptor 9; CMKAR3; IP10-R; CD183; MigR; GPR9; Interferon-Inducible Protein 10 Receptor; Chemokine (C-X-C Motif) Receptor 3; C-X-C Chemokine Receptor Type 3; IP-10 Receptor; Chemokine Receptor 3; CD183 Antigen; Mig Receptor; CXC-R3; CXCR-3; CD182; Mig-R

Background

Gene Name: CXCR3

NCBI Gene Entry: [2833](#)

UniProt Entry: [P49682](#)

Application Information

Clonality: Armenian Hamster monoclonal antibody

Clone ID: 24GB15680

Species Reactivity: Mouse

Applications Tested: Flow Cytometry (FCM)

Conjugation: FITC-Cy5

Immunogen

The immunogen information is proprietary.

Isotype

Armenian Hamster IgG

Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

Storage

Store at 4 °C and protected from prolonged exposure to light.

Recommended Dilutions

5 µL/Tests

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041


ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-CD183 Armenian Hamster Monoclonal Antibody, FITC-Cy5 Conjugated

PAGE 2



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM